

Tele ECG Solution

www.artemhealthtech.com



12 LEAD
SIMULTANEOUS



EASY TO
USE



COST
EFFECTIVE



TRUSTED BY
SPECIALIST



CLOUD BASED
ANY TIME/ANY WHERE



QUICK & 24 X 7
REPORTING



ALERTS FOR
CRITICAL CASES



5 LACS+
ECG REPORTED

“Artem HealthTech” 419, 4th floor, Block-D, Times Square Grand, Sindhu Bhavan Road, Thaltej, Ahmedabad-380 059, Gujarat, India.



+91 635 969 3050



info@artemhealthtech.com

Online ECG Reporting



ARTEM

HEALTH TECH

Formerly Known as Express Diagnostics Pvt. Ltd.

Care Provider: "Express ECG"



Name of Patient

Location of Patient/Hospital

Gender/DOB (Age) : Gender / DD-MM-YY (Y / M)


Medico ID : 14 Digit Unique ID

Referred By

Date : DD-MM-YY /HH-MM PM

History

REPORT ON ECG



Spring - 25mm/Sec - 10mm/mV. Date - Medical Hospital - V2.0.0.06



VITALS	TEMP : - (F)	PULSE RATE : - /MIN	RBS : - mg/dL
	HR : 128 /MIN	BP : 0 / 0 mmHg	SPO2 : 0.0 %
MEASUREMENTS*	PR : 145.0 ms	QT : 364.58 ms	P : 40.52 deg
(ECG Parameters)	ST : 0.48 ms	QTc : 531.33 ms	QRs : 66.17 deg
	R-R : 470.83 ms	QRS : 125.0 ms	T : 0.0 deg
FINDINGS	: SINUS TACHYCARDIA BROAD P WAVE NO SIGNIFICANT ST AND T WAVE CHANGES		
IMPRESSION	: SINUS TACHYCARDIA CONSIDER ATRIAL ENLARGEMENT		
RECOMMENDATION	: CLINICAL CORRELATION		

This is electronically authenticated report, hence doesn't require signature.

* Software calculated values; to be verified manually.
 Printed By : Name of Technician DD-MMM-YYYY / HH-MM-PM

Reported By
 Artem HealthTech HQ

Name of Reporting Doctor Reg. No. #####

How it Works?

- 1 Download the Express ECG Application from the Play Store.
- 2 Register the patient and apply the 12-lead ECG Belt / patient cable.
- 3 Connect the device, capture the ECG, & then upload it on ECG application.
- 4 Download the reported ECG or share it via WhatsApp, email, or print. In case of a critical report, connect with the specialist.



Product Specification

Make & Model	ET - U12
Channels	12 Channel Simultaneous
Cable Type	10 Lead ECG Chest Belt / Electrons & Chest Bulbs
Electrode Leads	Chest Belt with DB Connection
	0.05 to 100hz
Speed	25 mm/s 12.5
A/D Resolution	24 bits
	10 mm / mV
Sample Rate	500 Samples / Second Per Channel
Data Capturing Length	10 Sec. / 30 Sec. / 60 Sec.
Computing Device Required	PC / Laptop / Android Mobile & Tablet
Display	On Computing Device
Memory	Dependent On Computing Device (Automatic Data Store)
Software Application	Android App / Web App
Interface	USB / BLUE TOOTH
OTG Support	Yes
ECG Report	PDF Format A4 Normal Page
Dimension (L x W x H)	75mm x 55mm x 25mm



Artem's Tele ECG Salient Features

- 7 Low cost, and maintenance free
- 7 Compact and portable device
- 7 Powerful measurement & analysis function.
- 7 Compatible with any computing device.
- 7 Compliant to ABDM (all milestones)
- 7 Mobile app for technicians, reporting Doctors and Patient with EHR / EMR.
- 7 Multiple options to share and print reports
- 7 Facility to compare with previous ECG
- 7 Easy to use and easy to train solution
- 7 24 x 7 Online Reporting Service



REDEFINING CARE